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References:

1. Shane E, Burr D, Abrahamsen B, et al. Atypical subtrochanteric and diaphyseal femoral fractures: second report of a task force of the American Society for Bone and Mineral Research. *J Bone Miner Res.* 2014;29(1):1-23.
2. Brauer CA, Coca-Perraillon M, Cutler DM, Rosen AB. Incidence and mortality of hip fractures in the United States. *JAMA.* 2009;302(14):1573-1579.
3. Odvina CV, Zerwekh JE, Rao DS, Maalouf N, Gottschalk FA, Pak CY. Severely suppressed bone turnover: a potential complication of alendronate therapy. *J Clin Endocrinol Metab.* 2005;90(3):1294-1301.
4. Jha S, Wang Z, Laucis N, Bhattacharyya T. Trends in Media Reports, Oral Bisphosphonate Prescriptions, and Hip Fractures 1996–2012: An Ecological Analysis. *J Bone Miner Res.* 2015;30(12):2179-2187.
5. ASBMR, et al. *Call to Action to Address the Crisis in the Treatment of Osteoporosis.* September 19, 2016.
6. Kolata G. Fearing Drugs' Rare Side Effects, Millions Take Their Chances With Osteoporosis. *The New York Times.* June 1, 2016.
7. Kim SC, Kim DH, Mogun H, et al. Impact of the U.S. Food and Drug Administration's Safety-Related Announcements on the Use of Bisphosphonates After Hip Fracture. *J Bone Miner Res.* 2016;31(8):1536-1540.
8. Lewiecki EM. Hip Fractures and Declining DXA Testing: At a Breaking Point? *J Bone Miner Res.* 2016;31(Suppl 1).
9. Kim S, Yang KH, Lim H, et al. Detection of Prefracture Hip Lesions in Atypical Subtrochanteric Fracture with Dual-Energy X-ray Absorptiometry Images. *Radiology.* 2014;270(2):487-495.
10. Harborne K, Hazlehurst JM, Shanmugaratnam H, et al. Compliance with established guidelines for the radiological reporting of atypical femoral fractures. *Br J Radiol.* 2016;89(1057):20150443.

© 2021 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, enCORE, Prodigy and Lunar iDXA are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

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